

***InVitroGRO*[™] CP Rodent Medium**

Product No.	Description	Size
Z990028	<i>InVitroGRO</i> [™] CP Rodent Medium	250 mL
Z990027	<i>Torpedo</i> [™] Rodent Antibiotic Mix	5.5 mL

Product Description:

InVitroGRO CP Rodent Medium is used for plating BioreclamationIVT rat and mouse cryoplateable hepatocytes. *InVitroGRO* CP Rodent Medium is formulated without antibiotics to increase its shelf life. BioreclamationIVT recommends the use of the *Torpedo*[™] Rodent Antibiotic Mix (Z990027) with *InVitroGRO* CP Rodent Medium (Z990028).

Stability: 3 months for *InVitroGRO* CP Rodent Medium without *Torpedo* Rodent Antibiotic Mix
7 days for complete *InVitroGRO* CP Rodent Medium

Storage: 4°C

Caution: Treat all products containing human and monkey-derived materials a potentially infectious, as no known test methods can offer assurance that products derived from human or monkey tissues will not transmit infectious agents.

All products are for research use only. Do not use in animals or humans. These products have not been approved for any diagnostic or clinical procedures.