

InVitroGRO[™] HI Medium

Product No.	Description	Size
Z99009	<i>InVitroGRO</i> HI Medium	250 mL
Z990012	<i>InVitroGRO</i> HI Medium	500 mL
Z990002	<i>InVitroGRO</i> HI Medium	1 L

Product Description:

InVitroGRO HI Medium is a hepatocyte incubation medium used for culturing primary hepatocytes¹. *InVitroGRO* HI Medium is formulated without antibiotics to increase its shelf life. BioreclamationIVT recommends the use of the *Torpedo*[™] Antibiotic Mix with *InVitroGRO* HI Medium.

Stability: 6 months for *InVitroGRO* HI Medium without *Torpedo* Antibiotic Mix
7 days for Complete *InVitroGRO* HI Medium

Storage: 4°C

Related Products:

Product No.	Size
Z99000	<i>Torpedo</i> Antibiotic Mix 5.5 ml
Z990007	<i>Torpedo</i> Antibiotic Mix 11 mL
Z990008	<i>Torpedo</i> Antibiotic Mix 22 mL

Reference:

1. Roymans, D.; Van Looveren, C.; Leone, A.; Parker, J. B.; McMillan, M.; Johnson, M. D.; Koganti, A.; Gilissen, R.; Silber, P.; Mannens, G.; Meuldermans, W. Determination of cytochrome P450 1A2 and cytochrome P450 3A4 induction in cryopreserved human hepatocytes. *Biochem. Pharmacol.* **2004**, *67*(3), 427–437.

Caution: Treat all products containing human and monkey-derived materials a potentially infectious, as no known test methods can offer assurance that products derived from human or monkey tissues will not transmit infectious agents.

All products are for research use only. Do not use in animals or humans. These products have not been approved for any diagnostic or clinical procedures.