

## ***InVitroGRO*<sup>™</sup> KHB**

Product No.	Description	Size
Z99074	<i>InVitroGRO</i> KHB	250 mL

### **Product Description:**

*InVitroGRO* KHB is a modified Krebs-Henseleit Buffer for use with hepatocytes. BioreclamationIVT recommends the use of the *InVitroGRO* KHB along with our *InVitroGRO* hepatocyte media when using LiverPool<sup>™</sup> 5-, 10-, 20- and 50-Donor Pooled Cryopreserved Human Hepatocytes and regular Cryopreserved Human and Animal Hepatocytes.

**Stability:** 6 months  
**Storage:** 4°C

**Caution:** Treat all products containing human and monkey-derived materials a potentially infectious, as no known test methods can offer assurance that products derived from human or monkey tissues will not transmit infectious agents.

All products are for research use only. Do not use in animals or humans. These products have not been approved for any diagnostic or clinical procedures.